

# HIGHLIGHTS



# GC EUROPE



SMILES AROUND THE WORLD

**GC**



## GC G-ænial Universal Flo

### High-Viscosity Light-Cured Flowable Hybrid Composite

GC G-ænial Universal Flo presents a number of intelligent and convenient solutions that take it further than the standard flowables on the market: the physical performance of a regular composite material but in a flowable consistency. With GC G-ænial Universal Flo you get an exclusive combination of a flowing viscosity, for easy and convenient placement, combined with outstanding physical properties for safe and durable restorations.



## Characteristics & benefits

- Carefully balanced viscosity: smoothly flowing for easy placement but also highly thixotropic to hold its shape
- Very high gloss rate and gloss retention rate
- Can be considered to be a self polishing material
- Superb invisible aesthetics without any compromises
- Broad range of 15 shades and 3 different levels of translucencies
- Specially designed syringe that offers an ergonomic solution for improved application

## Packages

GC G-ænial Universal Flo (2 ml/3.4 g per syringe)	
004213	BW, 1 syringe
004202	A1, 1 syringe
004203	A2, 1 syringe
004204	A3, 1 syringe
004205	A3.5, 1 syringe
004206	A4, 1 syringe
004207	B1, 1 syringe
004208	B2, 1 syringe
004209	B3, 1 syringe
004210	C3, 1 syringe
004214	CV, 1 syringe
004211	AO2, 1 syringe
004212	AO3, 1 syringe
004215	JE, 1 syringe
004216	AE, 1 syringe

## GC G-ænial Bond

### One Component Self-Etching Light-Cured Adhesive

With the selective etching approach, GC G-ænial Bond offers the best of both worlds: the simplicity and reduced post-operative sensitivity of a self-etch adhesive together with the greater bond strength that was traditionally found only with etch & rinse adhesives. The dimethacrylate monomer increases its permeability into enamel and dentine, while the increased level of phosphate ester monomer optimizes etching. The resulting nano-interaction zone – and the HEMA-free formulation are key factors for the superior and long-lasting bond strength.



## Characteristics & benefits

- Can be used as self-etch adhesive or with selective etching of enamel
- Creation of a nano-interaction zone for superior bond strength to both enamel and dentine
- Excellent marginal integrity regardless of the technique and the composite material
- HEMA-free for a long-lasting bond
- Clinical reliability proven in independent tests

## Packages

GC G-ænial Bond	
004217	Starter Kit (1 x 5 ml bottle, 50 disposable applicators (fine), 20 x disposable dispensing dish)
004219	3-bottle Pack (3 x 5 ml bottle)
004220	Refill (1 x 5 ml bottle)

## GC MI Paste Plus

Remineralising protective crème with triple the benefit:  
Strengthen. Protect. Replenish.

MI Paste Plus is the ideal way to give your patients maximum all-round protection. Thanks to RECALDENT™\* – its revolutionary ingredient! MI Paste Plus remineralises enamel lesions due to the high level of calcium, phosphate and fluoride it provides. MI Paste Plus with CPP-ACP(F) is the only product that delivers the ideal calcium, phosphate and fluoride ratio of 5:3:1 and releases all three of the ions needed to form acid-resistant fluorapatite by BOTH remineralisation and fluoridation.



Recaldent



## Characteristics & benefits

- Buffers the pH changes in plaque
- Impairs the adhesion and growth of Streptococcus mutans and sobrinus
- Remineralises early enamel lesions and makes enamel more resistant
- Optimises the transport and intake of fluoride by enamel

## Indications

- During and/or after orthodontic treatment, especially on white spots
- For patients with medium or high susceptibility to caries and who need support due to very poor plaque control
- For patients who suffer from an acidic oral environment, erosion and gastric reflux
- For pregnant women

## Packages

GC MI Paste Plus (40 g, 35 ml per tube)	
002614	Assorted Pack 10 pcs, 2 of each flavour (Melon, Strawberry, Tutti Frutti, Mint and Vanilla)
002886	Strawberry, Pack of 10
002887	Melon, Pack of 10
002888	Vanilla, Pack of 10
002889	Tutti Frutti, Pack of 10
002621	Mint, Pack of 10

\* CPP-ACP was developed at the School of Dental Science at the University of Melbourne Victoria/Australia. RECALDENT™ is a trade mark used under license. RECALDENT™ (CPP-ACP) is derived from milk casein. Do not use on patients with a milk protein and/or hydroxybenzoates allergy.

## GC EQUIA Fil & Coat

### Self-Adhesive Posterior Restorative

The name GC EQUIA stands for “Easy – Quick – Unique – Intelligent – Aesthetic” and denotes a totally new approach to filling therapy: restorations based on glass ionomer technology have never been so aesthetic and translucent, high-performing and economical! GC EQUIA is the first to combine quick and easy handling with perfect physical and incomparably aesthetic properties. Recent clinical studies and evidences indicate that GC EQUIA can be used as a long term restorative material in stress bearing Class I and in limited sized Class II cavities! This makes GC EQUIA the most ideal alternative to amalgams and compomers.



## Characteristics & benefits

- No bonding agent needed, chemical adhesion with tooth structure
- Low moisture sensitivity
- Bulk placement with only 3'30" in total needed from start to finish
- Filler content provides wear resistance and fracture toughness
- Optimal protection of marginal seal for long-lasting restorations
- Tooth coloured restorative material, with real translucency and natural gloss and smoothness

## Packages

<b>GC EQUIA (mixed volume per capsule 0.14 ml)</b>			
<b>Intro Pack (box: 50 capsules EQUIA Fil, 4 ml EQUIA Coat &amp; accessories)</b>			
900603	Assorted	900558	B1
900554	A1	900559	B2
900555	A2	900560	B3
900556	A3	900561	C4
900557	A3.5		
<b>Promo Pack (2 boxes: 50 capsules EQUIA Fil, 4 ml EQUIA Coat &amp; accessories)</b>			
900562	Assorted & A3	900567	A3.5 & A3
900563	A1 & A3	900568	B1 & A3
900564	A2 & A2	900569	B2 & A3
900565	A2 & A3	900570	B3 & A3
900566	A3 & A3	900571	C4 & A3
<b>Refill (box with 50 capsules)</b>			
004259	Assorted	004264	B1
004260	A1	004265	B2
004261	A2	004266	B3
004262	A3	004267	C4
004263	A3.5		
<b>Refill EQUIA Coat</b>			
004268	Bottle of 4 ml, 20 Disposable Dispensing Dishes, 50 Micro-Tips Applicators, 1 Applicator Holder		

## GC EXA'lence

### Vinyl Polyether Silicone Impression Material

With the introduction of GC EXA'lence GC has been able to develop the ideal solution for all of your clinical situations and to remove the dilemmas involved in the choice of your impression material. The combination of intrinsic hydrophilicity and tear strength into one impression material ensures that from now on you will be able to take excellent impressions under every condition. The EXA'lence range sees GC meeting the demands of both dentists and technicians by providing a complete system with different viscosities and setting times. Achieving perfect impressions for perfectly fitting restorations under unpredictable circumstances is now possible for every dental professional.



2010  
REALITY  
★★★★

THE DENTAL  
ADVISOR  
★★★★

## Characteristics & benefits

- Intrinsic hydrophilicity ensures wettability of the material before setting for a perfectly detailed capture, even under wet conditions
- Remains hydrophilic after setting, for the highest precision during pouring of the model
- High tear strength eliminates the risk of damage during removal of the impression from the mouth and during removal of the model from the impression
- Long working time and short setting time
- Pleasing mild mint taste for the patient
- Easy to read colour contrast between the different viscosities

## Packages

GC EXA'lence	
Refill	Fast or Regular Set: 2 x 48 ml cartridges, Fast Set, 6 mixing tips: Extra Light Body, Light Body, Medium Body, Heavy Body, Heavy Body Rigid
GC EXA'lence Putty	
137810	Standard 1-1 Pack: 500 g (278 ml) base, 500 g (278 ml) catalyst & accessories
137866	Clinic 5-5 Pack: 5 x 500 g (278 ml) base, 5 x 500 g (278 ml) catalyst, accessories
GC EXA'lence 370	
Standard Kit	Fast or Regular Set: 370 ml cartridge, 10 dynamic mixing tips & 1 bayonet locking ring: Heavy Body or Monophase
Refill	Fast or Regular Set: 2 x 370 ml cartridges, Heavy Body or Monophase
Clinic Pack	Fast or Regular Set: 6 x 370 ml cartridges, Heavy Body or Monophase
Accessories	
150004	50 x dynamic mixing tips (yellow) for 370 ml
150005	1 x bayonet locking ring for 370 ml
138900	7 ml universal VPS adhesive bottle

## GC Initial IQ – One Body – Layering Concept

One Body, Layering-over-Metal and Layering-over-Zircon

Through the development of the GC Initial IQ – One Body – Layering Concept, GC can now offer a compact ceramic set for the fast and easy reproduction of the standard Vita® shades\*. This new system enables technicians to concentrate on the build-up and create vitality and a natural appearance with the GC Initial IQ Lustre Pastes NF in the final step. The GC Initial IQ – One Body approach stands for consistent quality and aesthetical results time and again with fewer materials.



## Characteristics & benefits

- Layering-over-Metal is compatible with standard porcelain alloys
- Layering-over-Zircon is compatible with all Zr-based technologies
- Crown build-up with only one powder allows full attention to shaping and function
- Fast reproduction of standard Vita® colours\* following the 'Paint by numbers' with GC Lustre Pastes NF
- No extra investments required
- Can be used as an excellent starter module for the Initial ceramic line

## Packages

<b>GC Initial IQ One Body, Layering-over-Metal (L-O-M)</b>			
877077	Set: 5 x MC/Paste Opaque 3 g (A2, A3, B2, C2, D2), 4 x IQ/ L-O-M Body 50 g (A, B, C, D), modelling liquid, Glass Plate, Brush, Technical Manual		
<b>GC Initial IQ One Body, Layering-over-Zircon (L-O-Z)</b>			
877079	Set: 7 x Light Reflective Liners 3 g (1-6, BL), Liner Thinner, 5x IQ/ L-O-Z Body 50 g (A, B, C, D, BL), modelling liquid, Glass Plate, Brush, Technical Manual)		
877080	Initial IQ, One Body, Light Reflective Liner Set (7 x Light Reflective Liners 3 g (1-6, BL), Liner Thinner)		
<b>GC Initial IQ, Lustre Pastes NF</b>			
877078	Set: 1 x Neutral, 1 x Value, 4 x Body Shades, 8 x Enamel effect shades, Diluting Liquid, Refresh Liquid, 1 mixing dish, 2 brushes, Picture IFU)		
<b>Refills</b>			
876401	Neutral, 4 g	876409	Enamel Effect Dark Grey, 4 g
876402	Body A, 4 g	876410	Enamel Effect Light Blue, 4 g
876403	Body B, 4 g	876411	Enamel Effect Dark Blue, 4 g
876404	Body C, 4 g	876412	Enamel Effect INcizio, 4 g
876405	Body D, 4 g	876413	Enamel Effect Olive, 4 g
876406	Enamel Effect Vanilla, 4 g	876414	Enamel Effect Value, 4 g
876407	Enamel Effect White, 4 g	876220	Diluting Liquid, 8 ml
876408	Enamel Effect Light Grey, 4 g	876400	Refresh Liquid, 8 ml
<b>Accessories</b>			
877051	Brush 00, 1 pce	877054	Plastic cover, 1 pce
877052	Brush 2, 1 pce	877055	Aluminium box, 1 pce
877053	Mixing dish, 1 pce		

\* Vita® is a registered trademark of Vita® - Zahnfabrik Bad Säckingen, Germany.



**GC EUROPE N.V.**

Head Office  
Researchpark  
Haasrode-Leuven 1240  
Interleuvenlaan 33  
B - 3001 Leuven  
Tel. +32.16.74.10.00  
Fax. +32.16.40.48.32  
info@gceurope.com  
<http://www.gceurope.com>

**GC BENELUX B.V.**

Edisonbaan 12  
NL - 3439 MN Nieuwegein  
Tel. +31.30.630.85.00  
Fax. +31.30.605.59.86  
info@benelux.gceurope.com  
<http://benelux.gceurope.com>

**GC FRANCE s.a.s.**

9 bis, Avenue du Bouton d'Or -  
BP 166  
F - 94384 Bonneuil sur Marne  
Cedex  
Tel. +33.1.49.80.37.91  
Fax. +33.1.49.80.37.90  
info@france.gceurope.com  
<http://france.gceurope.com>

**GC UNITED KINGDOM Ltd.**

12-15, Coopers Court  
Newport Pagnell  
UK - Bucks. MK16 8JS  
Tel. +44.1908.218.999  
Fax. +44.1908.218.900  
info@uk.gceurope.com  
<http://uk.gceurope.com>

**GC ITALIA S.r.l.**

Via Calabria 1  
I - 20098 San Giuliano Milanese  
Tel. +39.02.98.28.20.68  
Fax. +39.02.98.28.21.00  
info@italy.gceurope.com  
<http://italy.gceurope.com>

**GC IBÉRICA**

Dental Products, S.L.  
Edificio Codesa 2  
Playa de las Americas, 2, 1°, Of. 4  
ES - 28230 Las Rozas, Madrid  
Tel. +34.916.364.340  
Fax. +34.916.364.341  
info@spain.gceurope.com  
<http://spain.gceurope.com>

**GC AUSTRIA GmbH**

Swiss Office  
Bergstrasse 31  
CH - 8890 Flums  
Tel. +41.81.734.02.70  
Fax. +41.81.734.02.71  
info@switzerland.gceurope.com  
<http://switzerland.gceurope.com>

**GC AUSTRIA GmbH**

Tallak 124  
A - 8103 Rein bei Graz  
Tel. +43.3124.54020  
Fax. +43.3124.54020.40  
info@austria.gceurope.com  
<http://austria.gceurope.com>

**GC GERMANY GmbH**

Seifgrundstrasse 2  
D - 61348 Bad Homburg  
Tel. +49.61.72.99.59.60  
Fax. +49.61.72.99.59.66.6  
info@germany.gceurope.com  
<http://germany.gceurope.com>

**GC NORDIC AB**

Danish Branch  
Dampfærgevej 3  
DK - 2100 København  
Tel. +45 23 26 03 82  
info@denmark.gceurope.com  
<http://nordic.gceurope.com>

**GC NORDIC AB**

Finnish Branch  
Vanha Hommaksentie 11B  
FIN - 02430 Masala  
Tel. & Fax. +358.9.221.82.59  
info@finland.gceurope.com  
<http://finland.gceurope.com>

**GC EUROPE N.V.**

East European Office  
Siget 19b  
HR - 10020 Zagreb  
Tel. +385.1.46.78.474  
Fax. +385.1.46.78.473  
info@eeo.gceurope.com  
<http://eeo.gceurope.com>